



**VISIONEC-LAB KFT.**

Safety Caps and much more...

## Certificate of Compliance

### Vision Shut-off valve

#### Applies to Part Numbers:

VLV-SF1, VLV-SF1-03, VLV-SF1-10

During a 40-hour working week, under regular laboratory conditions (20%<RH<70%, T=20-25°C, Flowrate<2ml/min) **Vision Eluent Shut-off Valves** can be used effectively to prevent solvent vapors escaping from the storage apparatus and to filter the incoming air for more than 6 months but due to the variety of other factors like: different solvent type/composition, different volatile eluent modifiers, pressure, temperature, humidity, dust in the laboratory air we do recommend to replace the filters after the recommended 6-month period to provide flawless operation.

**Lot Number:** \_\_\_\_\_

**Leak test:** \_\_\_\_\_

**Initial:** \_\_\_\_\_

**Date of issue:** \_\_\_\_\_



info@visionec.hu



2132 Göd, Attila utca 13.



www.visionsafetycap.com - www.visionhpfittings.com